


**EXPERIMENT NO:**

**DATE:**

**PRINCIPLE:** Hardness of water is mainly due to the presence of calcium and magnesium salts in it. Total hardness is the sum of temporary hardness (due to bicarbonates of calcium and magnesium) and permanent hardness (due to chlorides, sulphates etc, of calcium and magnesium). Ethylenediaminetetraacetic acid is a reagent, which reacts with metal ions like  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$  forming complex compounds. Therefore, this reagent can be used to determine the concentration of hardness causing substances.



The completion of the reaction (end point of the titration) is identified using Erichrome black T indicator. This is an organic dye, blue in colour. It also forms relatively less stable complexes with bivalent metal ions of Ca, Mg etc., which produces wine-red colour. Therefore, addition of the indicator to hard water produces wine-red colour. When EDTA is added to hard water, it first reacts with free metal ions and then attacks the metal-indicator complex. The latter reaction can be represented as



So, at the end point a change from wine-red to blue colour is observed. Since the reaction involves the liberation of  $H^+$  ions and the indicator is sensitive to the concentration of  $H^+$  ions (metal indicator complex is stable at pH 10) of the solution, a constant pH of around 10 has to be maintained. For this purpose ammonia-ammonium chloride buffer is used.

### Part-A:Preparation of standard EDTAsolution:

Weigh the weighing bottle containing disodium salt into the funnel placed on a 250cm<sup>3</sup> volumetric flask. Weigh the bottle again. The difference between the two weights will give the amount of EDTA transferred. Pour small quantities of water over the salt on the funnel and transfer the salt into the flask. Wash the funnel with the same water 3-4 times. Dissolve the salt by adding about 100-150 cm<sup>3</sup> of deionised water. Make up the solution to the mark and shake well for uniform concentration.

Pipette out 25cm<sup>3</sup> of the given sample of hard water into a clean conical flask. Add 3cm<sup>3</sup> of NH<sub>3</sub>-NH<sub>4</sub>CL buffer followed by 5-6 drops of Erichrome black-T indicator. Titrate this against EDTA taken in a burette till the colour changes from wine-red to blue. Note down the burette readings and repeat the titration to get concordant values.

Experiment no

## Observation and Calculations

### Part A: Preparation of standard EDTA solution

Weight of weighing bottle and EDTA salt:..... g

Weight of empty weighing bottle:.....g

Weight of EDTA salt:.....g

$$\text{Molarity of EDTA: } \frac{\text{Weight of EDTA salt} \times 4}{\text{Molecular weight of EDTA}} = \frac{W \times 4}{372.24} \quad (X)$$

### Part B: Determination of total hardness of water sample

Burette: Standard EDTA solution

Conical flask: 25 cm<sup>3</sup> of water sample + 3 cm<sup>3</sup> of NH<sub>4</sub>Cl-NH<sub>4</sub>OH buffer solution

Indicator: A pinch of Eriochrome Black T indicator

End point: Wine red to Blue

Trial	I	II	III
Burette readings			
Final reading			
Initial reading			
Volume of EDTA run down in cm <sup>3</sup>			

1000cm<sup>3</sup> of 1M EDTA solution = 100g of CaCO<sub>3</sub>

Therefore Vcm<sup>3</sup> of X M EDTA =  $\frac{X \times V \times 100}{1000}$

= (a g) of CaCO<sub>3</sub>

25cm<sup>3</sup> of hard water contains ..... g of CaCO<sub>3</sub>

Therefore total hardness of water =  $\frac{a \times 10^6}{25}$  = ..... =ppm

**RESULT: Total hardness of water=**

**ppm**